



We care for your back health



# Low back pain is identified as the leading cause of disability in the world.

Low back pain is the leading cause of disability and loss of quality of life worldwide, according to a study published in the journal *Annals of the Rheumatic Diseases*.

**9,4%** According to the study's data, **9,4 % of the world population suffers from low back pain**, with the highest prevalence in Western European countries.

**Ageing** The authors warn: **“given the ageing population throughout the world, and especially in developing countries, the number of people with low back pain will increase substantially in the coming decades”**.

**Lower quality of life** Of a list of 290 diseases included in the study, **low back pain ranks first** in terms of lost years of disability.



## Men, the most at risk

The study identifies men as the people with the highest risk of suffering from low back pain, aged between 35 and 65 years. In the case of farmers, the risk of low back pain is up to 3.7 times higher than that of the rest of the workers.



**Primpine** Allies with your well-being



# Primspine *Allies with your well-being*

We know that each person has different needs. That's why **Primspine** has created a wide range of solutions for your back, using exclusive technologies and fabrics that will act effectively on your pain or injuries.

Discover our new line of orthoses  
*offering specific solutions for your needs.*





# Exclusive fabrics.

## Technology that cares for your back

### The fitting: the secret to comfort

**A key feature with a direct effect on the level of stabilization of the orthosis that multiplies the feeling of comfort.**

**PRIM** always at the forefront, incorporates in its **PrimSpine** range an exclusive elastic fabric, with a double channel structure, highly breathable, which enhances the adaptation to the body, providing greater compression and multiplying the sensation of comfort.

Exclusive double-channel boostfit fabric

## BOOST fit

+ Fitting  
+ Compression  
+ Comfort



Scan me and discover our exclusive fabric



## fitCONTOUR

Inelastic three-layer fitContour fabric

### Capable of adapting to all anatomies.

The openings in the upper layer produce a "chimney effect" that allows the skin to breath.

**Primspine** Allies with your well-being



new  
line

# Primspine MOVE+

*Allies with your well-being*

**Contention & protection**



PRS607/PRS606



## Lumbar support belt

Made of elastic monofilament fabric with 2 rear stays that provide light support at the lumbar level. Incorporates a Velcro®-type front closure to allow adjustments to the abdominal compression.

# MOVE+




### Indications:

mild low back pain requiring light support, postural rehabilitation, muscular hypotonia, small and medium hernias.

PRS606 ●

Size/Size range	A cm.	B cm.
S 75-95	24	24
M 95-115	24	24
L 115-130	24	24
XL 130-140	24	24



Measurement: abdomen circumference

A: Front height  
B: Back height





*Allies with your well-being*

# Primospine CARE

**Stabilization level 1**



PRS610



## Semi-rigid lumbosacral back support belt Lumbar criss-cross

Made with **boostFit** fabric, it features a lumbar criss-cross and 4 rear stays to provide a high level of support. The front microvelcro closure allows for a secure, customised level of abdominal compression.

# CARE



**BOOST fit**

**1** Stabilization level

### Indications:

mild lower back pain, muscle weakness, mild lumbosciatica.

PRS610BG ● PRS610G ●

High-frequency technology

Breathable Latex-free

Size/Size range	A cm.	B cm.	A cm.	B cm.
XS 70-80	17	26	17	26
S 77-90	17	26	17	26
M 87-100	17	26	17	26
L 95-110	17	26	17	26
XL 107-120	17	26	17	26
XXL 117-130	17	26	17	26



A: Front height  
B: Back height

Measurements: Circumference at height of pelvis





*Allies with your well-being*

# Primpine CARE+

*Stabilization level II*



PRS620

## Semi-rigid lumbosacral back support belt Double lumbar criss-cross tensioner

Made with **boostFit** fabric, it features 4 lumbar stays and a double criss-cross tensioner at the back to strengthen and regulate compression to the area. The front microvelcro closure allows for a secure, customised level of abdominal compression.



**BOOST fit**

**2** Stabilization level

### Indications:

moderate lower back pain, moderate lumbosciatica, rehabilitation, muscle weakness, strains, sprains, early-stage degenerative lumbar spine processes, disc protrusion.



High-frequency technology



Transpirable Sin látex

### PRS620BG ● PRS620G ●

Size/Size range	A cm.	B cm.	A cm.	B cm.
XS 70-80	17	26	17	26
S 77-90	17	26	17	26
M 87-100	17	26	17	26
L 95-110	17	26	17	26
XL 107-120	17	26	17	26
XXL 117-130	17	26	17	26



A: Front height  
B:

Measurements: Circumference at height of pelvis

*Allies with your well-being*



# Primospine

## MODULAR

**Stabilization level III**



PRS630/PRS631



## Semi-rigid lumbosacral back support belt Double tensioner band

Made with **boostFit** fabric, it features 4 malleable stays in an inelastic back support. Its double tensioner band, with double microvelcro closure at the front, provides both lumbar and abdominal compression. Incorporates a lumbar plate insert.

# MODULAR

Can be combined with the **dorsolumbar shoulder support** and the **PrimSpine Modular lumbosacral module**



**BOOST fit**

**3** Stabilization level

### Indications:


lower back pain, lumbosciatica, osteoarthritis, early-stage osteoporosis, osteomalacia, disc protrusion, postoperative recovery, muscle weakness, sprain.

**PRS630BG** ●  
**PRS630G** ●  
Baja

	Size/Size range	A cm.	B cm.
XS	70-80	17	26
S	77-90	17	26
M	87-100	17	26
L	95-110	17	26
XL	107-120	17	26
XXL	117-130	17	26

**PRS631BG** ●  
**PRS631G** ●  
Alta

	Size/Size range	A cm.	B cm.
XS	70-80	20	30
S	77-90	20	30
M	87-100	20	30
L	95-110	20	30
XL	107-120	20	30
XXL	117-130	20	30



Measurements: Circumference at height of pelvis

A: Front height  
B: Back height